ISOGV303

User Information

MICRO POWDERED VINYL GLOVES NON STERILE

Available size : XS/6, S/7, M/8, L/9, XL/10

Manufacturer Address :

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1 - EU Type-Examination

a - This product is classed as Category III Personal Protective Equipment (PPE) according to PPE Regulation (EU) 2016/425 and has been shown to comply with this Regulation through the Harmonised European Standards EN 420:2003+A1:2009, EN ISO 374-1:2016 and EN ISO 374-5:2016.

b - Notified Body responsible for certification and Module B compliance is SATRA Technology Europe Ltd(2777), Bracetown Business Park Clonee, Dublin D15 YN2P Ireland.

c - Notify Body responsible for internal production control plus supervised product checks at random intervals (Module C2) is SATRA Technology Europe Ltd(2777), Bracetown Business Park Clonee, Dublin D15 YN2P Ireland.

d - The EU Declaration of Conformity is accessible at www.mediprotec.com

2 - Marking

a - Micro Organism Hazards Pictogram : EN 374-2:2014 Resistance to Penetration. No leak detected during water leak test. Minimum AQL is 1.5 or EN performance level 2 ; Additional information obtainable from the manufacturer.


Performance Level	AQL	Inspection Levels
Level 3	<0.65	G1
Level 2	<1.5	G1
Level 1	<4.0	S4

b - Micro Organism Hazards Pictogram : EN 374-5:2016 Protection against Bacteria, Fungi and Virus.

Protection against bacteria and fungi : PASS

Protection against virus : PASS

The following pictogram is applied :


c - Chemical Hazards Pictogram : BS EN 16523-1:2015 Additional information on chemical resistance obtainable from manufacturer

EN ISO 374-1:2016 permeation levels are based on breakthrough times as follows :

Performance Level	1	2	3	4	5	6
Minimum break-through time (mins)	10	30	60	120	240	480

This product complies with Type C requirements and the following pictogram shall be used with reference to clause 6.3 of ISO 374-1.

3 - Performance and Limitation of Use

a - This product has been tested in accordance with EN 374-2:2014.

Water leak test : No leak detected

Air leak test : No leak detected

b - Gloves had been tested in accordance with BS EN 16523-1:2015 resistance to permeation by chemicals and achieved the following performance levels :

Chemicals	Performance Level
Diethylamine	The samples tested did not meet with the minimum breakthrough time for a performance level 1 to be achieved
40% Sodium Hydroxide (K)	6
96% Sulphuric acid (L)	1

WARNING

- This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.

- The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm - where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical used in a mixture.

- It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.

- when used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves;

- Before use, inspect the gloves for any defect or imperfections.

c - This product had been tested in accordance with BS EN 374-4:2013 and achieved the following degradation results :

Chemicals	Mean Degradation / %
Diethylamine	-45.9
40% Sodium Hydroxide (K)	-28.6
96% Sulphuric acid	-62.5

EN 374-4:2013 degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemicals.

d - This product provides protection against Bacteria, Fungi and Virus. The gloves had been tested in accordance with ISO 16604:2014 to meet the requirements of BS EN ISO 374-5:2016 for resistance to penetration by blood-borne pathogens-test method using Phi-X174 bacteriophage.

- The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

e - The gloves were found to meet with the REACH annex XVII requirements for Polycyclic Aromatic Hydrocarbons (PAHs).

f - Components used in glove manufacturing may cause allergic reactions in some users. If allergic reactions occur, seek for medical advice immediately.

4 - EN 420:2003+A1:2009

Clause 5.1 - Length and fit : PASS

Clause 5.2 - Dexterity : Level 5

Glove size	7 - S	8 - M	9 - L	10 - XL
Minimum length / mm	230	240	250	260

Performance level	1	2	3	4	5
Diameter of dexterity pin / mm	11.0	9.5	8.0	6.5	5.0

Gloves for Special Applications (EN 420:2003, Clause 5.1.3)

These gloves are manufactured with unique standard sleeve length .Therefore , they are not meeting the Clause 5.1.3 from rule EN420:2003 , describing the sleeves length for each size of gloves. Nevertheless, as special application for these gloves is to protect fingers and palm area from chemical splashes , they can be considered as totally adapted to their special use.

5 - Product Instruction for Use

a - Usage - For Single Use only. If re-used, the risk of contamination and infection increases due to improper cleaning processes; and increased risk of holes and tear during re-use to weakening of gloves by cleaning processes.

b - Sizing - Select the right size glove for your hand

c - Donning - Hold glove by the bead with one hand. Align the glove thumb with your other hand thumb and slide your hand into the glove, one finger into each glove finger. Pull by the glove palm to get a good fit. Don the other glove by the same procedure.

d - Inspection - Punctures or tears may occur after donning. Inspect each glove after donning, and immediately discontinue use if found damaged.

e - Doffing - Hold glove bead and pull toward the fingers until glove comes off.

f - Disposal - Properly disposal of all used gloves. Follow your Institution's policies for disposal.

6 - Handling and Storage

Store in a cool and dry place. Opened boxes should be kept away from fluorescent and sunlight. Gloves are packed in dispenser which is suitable for transport. Keep the gloves in the box when not in use.

7 - Shelf life

The shelf life of product is 5 years from date of manufacture.

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